



Product Specification
Micro sample tube 100 µl brilliant, 100 µl, screw cap

Page 1



Product description

| | |
|---------------------|--|
| Order number | 42.117 |
| Product description | Micro sample tube 100 µl brilliant, preparation: 1% brilliant cresyl blue solution, nominal volume: 100 µl, (LxØ) with cap: 46 x 10.8 mm, screw cap, cap: blue, conical base, 50 piece(s)/case |

Product characteristics

| | |
|---------------------------|-----------------------------------|
| Type of collection | Venous and capillary |
| Type of preparation | 1% brilliant cresyl blue solution |
| Preparation concentration | 1% brilliant cresyl blue solution |
| Closure type | screw cap |
| Example of use | Reticulocyte counting |
| Base shape | conical base |
| Volume of work | 100 µl |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Micro sample tube 100 µl brilliant, 100 µl, screw cap

Page 2

Size

| | |
|----------------------|---------|
| Sample volume | 100 µl |
| Diameter | 10.8 mm |
| Length of product | 46 mm |
| Length including cap | 46 mm |
| Length excluding cap | 44 mm |

Material & colours

| | |
|-------------------|--------------------|
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | blue |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 117000 |
| Depth of case | 283 mm |
| Width of case | 173 mm |
| Height of case | 135 mm |
| Case volume | 0.0066 cbm |
| Weight of product | 0.0015 kg |
| Weight of case | 0.92 kg |
| EAN of inner box | 4038917080897 |
| EAN case | 4038917365130 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

