



Product Specification
**Absorbent liner, suitable for Mailing container 126
x 30 mm, (LxW): 75 x 90 mm**

Page 1



Product description

Order number	78.670
Product description	Absorbent liner, suitable for Mailing container 126 x 30 mm, (LxW): 75 x 90 mm, 1,000 piece(s)/case

Product characteristics

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	Mailing container 126 x 30 mm

Size

Width of product	90 mm
Length of product	75 mm

Material & colours

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Absorbent liner, suitable for Mailing container 126
x 30 mm, (LxW): 75 x 90 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 4000
Type of smallest case
subpackaging
Piece(s) / inner box 1000
Piece (s) / outer case 4000
Piece(s) / pallet 72000
Depth of case 515 mm
Width of case 360 mm
Height of case 360 mm
Case volume 0.0667 cbm
Weight of product 0.001 kg
Weight of case 4.36 kg
EAN of inner box 4038917377492
EAN case 4038917159937

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

