




Product Specification
Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm

Page 1



Product description

Order number	78.671
Product description	Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm, 100 piece(s)/case

Product characteristics

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	cold transport container

Size

Width of product	55 mm
Length of product	80 mm

Material & colours

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm

Page 2

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2500
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	2500
Piece(s) / pallet	270000
Depth of case	383 mm
Width of case	193 mm
Height of case	190 mm
Case volume	0.014 cbm
Weight of product	0.0009 kg
Weight of case	2.57 kg
EAN of inner box	4038917083263
EAN case	4038917079860

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT