



## **Product description**

Order number	93.1650.010
Product description	Insertion plate, Version 1, material: PP, 1 piece(s)/case

## **Product characteristics**

Suitable for	Safety Tray (93.1650) & e.g. 1x Multi-Safe twin plus 500 (77.3884.005/.105), 1x rack with handle (93.851 - 93.851.1xx), 2x dispenser boxes (95.813.001), 3x multi-purpose cup (75.1355), 1x disinfection spray
--------------	--

## **Size**

Width of product	272 mm
Length of product	420 mm

## **Material & colours**

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1440
Case volume	0.0009 cbm
Weight of product	0.152 kg
Weight of case	0.31 kg
EAN of inner box	4038917132695
EAN case	4038917076715

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

