



Product description

Order number 93.1650.030
Product description Insertion plate, Version 3, material: PP, 1 piece(s)/case

Product characteristics

Suitable for Safety Tray (93.1650) & e.g. 1x Multi-Safe twin plus 500 (77.3884.005/.105), 2x dispenser boxes (95.813.001), 8x multi-purpose cup (75.1355), 1x disinfection spray

Size

Width of product 272 mm
Length of product 420 mm

Material & colours

Product material Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1440
Case volume	0.0009 cbm
Weight of product	0.172 kg
Weight of case	0.33 kg
EAN of inner box	4038917132718
EAN case	4038917076739

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

