



### **Product description**

Order number	95.813.001
Product description	Dispenser box, material: PS, lid assembled

### **Product characteristics**

Base shape	flat base
------------	-----------

### **Size**

Diameter of base	81 mm
Diameter of lid	83 mm
Diameter of opening	81 mm
Length of product	110 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	white
Closure material	Polystyrene (PS)
Colour of cap	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category                      no medical device | IVD  
Batched                                    no

Packaging

Minimum order qty.                    35  
Piece(s) / inner box                   35  
Piece (s) / outer case                 35  
Piece(s) / pallet                        1260  
Depth of case                           398 mm  
Width of case                            398 mm  
Height of case                           300 mm  
Case volume                             0.0475 cbm  
Weight of product                      0.0623 kg  
Weight of case                          2.77 kg  
EAN of inner box                        4038917121644  
EAN case                                 4038917074797

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

