



Product description

| | |
|---------------------|--|
| Order number | 95.963.007 |
| Product description | Wall holder, for Multi-Safe twin plus, material: stainless steel |

Product characteristics

| | |
|----------------|---------------------------------------|
| Example of use | disposal of pointed and sharp objects |
| Suitable for | Multi-Safe twin plus |

Size

| | |
|-------------------|--------|
| Width of product | 186 mm |
| Length of product | 165 mm |

Material & colours

| | |
|------------------|-----------------|
| Product material | stainless steel |
|------------------|-----------------|

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Latex-free | yes |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 20 |
| Piece(s) / pallet | 2880 |
| Depth of case | 323 mm |
| Width of case | 213 mm |
| Height of case | 145 mm |
| Case volume | 0.01 cbm |
| Weight of product | 0.246 kg |
| Weight of case | 5.24 kg |
| EAN of inner box | 4038917900997 |
| EAN case | 4038917901000 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

