



Product description

Order number	95.963
Product description	Wall holder, for Multi-Safe twin plus, euroMatic® , with adhesive/flat adapter assembled, material: stainless steel

Product characteristics

Design	adhesive/flat adapter assembled
Suitable for	Multi-Safe twin plus, euroMatic®

Size

Width of product	51 mm
Length of product	139 mm

Material & colours

Product material	stainless steel
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	30
Piece(s) / pallet	2880
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.1801 kg
Weight of case	5.81 kg
EAN of inner box	4038917084697
EAN case	4038917003469

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

