



Product description

| | |
|---------------------|---|
| Order number | 77.1515 |
| Product description | Disposal container, 2 l, transparent, label/print: yellow/black, material: PP, opening: round, height: 245 mm, flat cap, Ø of cap: 82 mm, cap: yellow, 30 piece(s)/carton |

Product characteristics

| | |
|-----------------------|--------------|
| Colour of print/label | yellow/black |
| Cap | flat cap |
| Closure type | screw cap |
| Shape of opening | round |

Size

| | |
|---------------------|--------|
| Max. volume | 2 l |
| Diameter | 116 mm |
| Diameter of base | 105 mm |
| Diameter of lid | 82 mm |
| Diameter of opening | 71 mm |
| Height of product | 245 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Colour of cap | yellow |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 30 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 30 |
| Piece (s) / outer case | 30 |
| Piece(s) / pallet | 360 |
| Depth of case | 608 mm |
| Width of case | 368 mm |
| Height of case | 515 mm |
| Case volume | 0.1152 cbm |
| Weight of product | 0.125 kg |
| Weight of case | 4.5 kg |
| EAN of inner box | 4038917086813 |
| EAN case | 4038917005319 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

