



## Product description

Order number	77.3884.005
Product description	Disposal container, Multi-Safe twin plus, 500 ml, yellow, round, with plastic label, label/print: yellow/red, material: PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 162 mm, Ø of cap: 89 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 126 piece(s)/carton
PZN number	16791920

## Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Max. volume	500 ml
Diameter	89 mm
Diameter of base	57 mm
Diameter of lid	89 mm
Height of product	162 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Closure material	Polypropylene (PP)
Colour of cap	red
Colour of lid	red

### Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order qty.	126
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	126
Piece(s) / pallet	3024
Depth of case	266 mm
Width of case	451 mm
Height of case	599 mm
Case volume	0.0719 cbm
Weight of product	0.0762 kg
Weight of case	10.36 kg
EAN of inner box	4038917542425
EAN case	4038917542418

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

