



Product description

Order number	77.3884.025
Product description	Disposal container, Multi-Safe twin plus, 2,500 ml, yellow, with plastic label, label/print: yellow/red, material: PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 197 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 70 piece(s)/carton

Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	2500 ml
Diameter of base	112 mm
Height of product	197 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Closure material	Polypropylene (PP)
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	70
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	70
Piece(s) / pallet	840
Depth of case	800 mm
Width of case	400 mm
Height of case	380 mm
Case volume	0.1216 cbm
Weight of product	0.222 kg
Weight of case	17.15 kg
EAN of inner box	4038917903288
EAN case	4038917903295

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

