



**Product description**

Order number	77.3884.130
Product description	Disposal container, Multi-Safe twin plus, 3,000 ml, yellow, with plastic label, label/print: yellow/red, biohazard labeling, material: PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 233 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 62 piece(s)/carton

**Product characteristics**

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Size**

Max. volume	3000 ml
Diameter of base	112 mm
Height of product	233 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

**Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	62
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	62
Piece(s) / pallet	744
Depth of case	800 mm
Width of case	400 mm
Height of case	370 mm
Case volume	0.1184 cbm
Weight of product	0.241 kg
Weight of case	16.53 kg
EAN of inner box	4038917903523
EAN case	4038917903530

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

