



## **Product description**

Order number	77.3898.111
Product description	Disposal container, Multi-Safe medi, 11 l, yellow, round, material: PP, with assembled bracket, opening: round, height: 220 mm, flat cap, Ø of cap: 290 mm, cap: red, TRBA 250, DIN EN ISO 23907, 250 piece(s)/carton
UN identification	no

## **Product characteristics**

Cap	flat cap
Design	with assembled bracket
Shape of opening	round
Removal aid - turning aid	Removal aid

## **Size**

Max. volume	11 l
Diameter of base	257 mm
Diameter of lid	290 mm
Diameter of opening	88 mm
Height of product	220 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

### **Purity & certification**

Satisfies the requirement	TRBA 250, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	250
Type of smallest subpackaging	pallet
Piece(s) / inner box	1
Piece (s) / outer case	250
Piece(s) / pallet	250
Depth of case	792 mm
Width of case	1090 mm
Height of case	1192 mm
Case volume	1.029 cbm
Weight of product	0.6135 kg
Weight of case	163.32 kg
EAN of inner box	4038917409513
EAN case	4038917409506

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

