



Product Specification

Urine container, 3 l, with inspection strip and enclosed label with instructions for use, brown, with light protection, graduated

Page 1



Product description

Order number	77.575.001
Product description	Urine container, Urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, brown, with light protection, graduated, material: PE, with paper label, english, cap: yellow, with inspection strip and enclosed label with instructions for use, 30 piece(s)/case

Product characteristics

Label/ Print	with paper label
Graduation	yes
Closure type	screw cap
Light protection	yes
Design	with inspection strip and enclosed label with instructions for use
Application	Urine collection
Viewing window	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Urine container, 3 l, with inspection strip and enclosed label with instructions for use, brown, with light protection, graduated

Page 2

Size

Sample volume	3.000.00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm
Length excluding cap	243 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	yellow

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	270
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm
Case volume	0.1335 cbm
Weight of product	0.1136 kg
Weight of case	4.21 kg
EAN of inner box	4038917266222
EAN case	4038917266215

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT