



Product Specification

**UriSet 24, without preparation, 3 l, with
Urine-Monovette® 10 ml, brown, with light
protection, graduated**

Page 1



Product description

Order number	77.578.010
Product description	UriSet 24, Complete set for optimum 24-hour urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, without preparation, brown, with light protection, graduated, material: PE, cap: yellow, with Urine-Monovette® 10 ml, 30 piece(s)/case

Product characteristics

Type of preparation	without preparation
Graduation	yes
Closure type	screw cap
Light protection	yes
Design	with Urine-Monovette® 10 ml
Application	Complete set for optimum 24-hour urine collection
Viewing window	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**UriSet 24, without preparation, 3 l, with
Urine-Monovette® 10 ml, brown, with light
protection, graduated**

Page 2

Size

Sample volume	3.000.00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm
Length excluding cap	243 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	yellow

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	270
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm
Case volume	0.1335 cbm
Weight of product	0.1846 kg
Weight of case	6.34 kg
EAN of inner box	4038917113861
EAN case	4038917080019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT