



**Product Specification**  
**Medicine cup, 30 ml, (LxØ): 40 x 37 mm,**  
**graduated, PP, transparent**

**Page 1**



**Product description**

Order number	75.1337.500
Product description	Medicine cup, max. working volume: 30 ml, (LxØ): 40 x 37 mm, Ø opening: 37 mm, transparent, graduated, material: PP, 75 piece(s)/bag

**Product characteristics**

Graduation	yes
Volume of work	30 ml

**Size**

Sample volume	30 ml
Diameter	37 mm
Diameter of opening	37 mm
Height of product	40 mm
Length of product	40 mm
Length excluding cap	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Medicine cup, 30 ml, (LxØ): 40 x 37 mm,**  
**graduated, PP, transparent**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	Medical device
CE certified	CE supplier
Batched	yes

**Packaging**

Minimum order qty.	3750
Type of smallest subpackaging	bag
Piece(s) / inner box	3750
Piece (s) / outer case	3750
Piece(s) / pallet	180000
Depth of case	190 mm
Width of case	390 mm
Height of case	390 mm
Case volume	0.0289 cbm
Weight of product	0.0015 kg
Weight of case	6 kg
EAN of inner box	4038917127721
EAN case	4042836001994

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-02

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**