



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 75.560.901   |
| Product description | Urine container, max. working volume: 148 ml, (ØxH): 73 x 67 mm, white, material: PS |

### **Product characteristics**

|                |        |
|----------------|--------|
| Graduation     | no     |
| Volume of work | 148 ml |

### **Size**

|                   |        |
|-------------------|--------|
| Sample volume     | 148 ml |
| Diameter          | 73 mm  |
| Height of product | 67 mm  |

### **Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | white            |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

**Packaging**

|                        |               |
|------------------------|---------------|
| Minimum order qty.     | 2500          |
| Piece(s) / inner box   | 100           |
| Piece (s) / outer case | 2500          |
| Piece(s) / pallet      | 90000         |
| Depth of case          | 457 mm        |
| Width of case          | 368 mm        |
| Height of case         | 356 mm        |
| Case volume            | 0.0599 cbm    |
| Weight of product      | 0.0028 kg     |
| Weight of case         | 8.16 kg       |
| EAN of inner box       | 4038917318198 |
| EAN case               | 4038917916547 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

