



Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown

Page 1



Product description

Order number	75.562.011
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, brown, with light protection, graduated, material: PP, screw cap, cap: brown, cap enclosed, 50 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Light protection	yes
Application	collection and storage of urine
Volume of work	100 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown

Page 2

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	7500
Depth of case	443 mm
Width of case	323 mm
Height of case	330 mm
Case volume	0.0472 cbm
Weight of product	0.0153 kg
Weight of case	4.54 kg
EAN of inner box	4038917113618
EAN case	4038917014038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

