




**Product Specification**  
**Snap-on lid, without spout, PE, red, 1,000**  
**piece(s)/bag**

**Page 1**



**Product description**

Order number	76.1414
Product description	Snap-on lid, without spout, suitable for 125-ml container 75.1355, material: PE, red, 1,000 piece(s)/bag

**Product characteristics**

Closure type	snap-on lid
Design	without spout
Suitable for	125-ml container 75.1355

**Material & colours**

Product material	Polyethylene (PE)
Colour of product	red

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Snap-on lid, without spout, PE, red, 1,000**  
**piece(s)/bag**

**Page 2**

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	390 mm
Width of case	265 mm
Height of case	590 mm
Case volume	0.061 cbm
Weight of product	0.0024 kg
Weight of case	3.14 kg
EAN of inner box	4038917086554
EAN case	4042836001499

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**