




**Product Specification**  
**Snap-on lid, for 100-ml container 75.570, PS,  
transparent**



**Product description**

Order number	76.571
Product description	Snap-on lid, for 100-ml container 75.570, material: PS, transparent, 100 piece(s)/tube bag

**Product characteristics**

Closure type	snap-on lid
Suitable for	100-ml container 75.570

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Snap-on lid, for 100-ml container 75.570, PS,**  
**transparent**

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	tube bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Case volume	0.0392 cbm
Weight of product	0.0009 kg
Weight of case	1.46 kg
EAN of inner box	4038917082570
EAN case	4038917001656

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

