



Product Specification

**Midstream urine kit, 70 ml, (LxØ): 115 x 95 mm,
for 70 ml container, graduated, PP, translucent,
with paper label**



Product description

Order number	51.9923.812
Product description	Midstream urine kit, max. working volume: 70 ml, (LxØ): 115 x 95 mm, Ø opening: 48.45 mm, for 70 ml container, translucent, graduated, material: PP, screw cap, cap enclosed, with paper label, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with paper label
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Application	collection and storage of urine
Example of use	Urine Collection through funnel into 70ml container
Suitable for	70 ml container
Base shape	flat base
Volume of work	70 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Midstream urine kit, 70 ml, (LxØ): 115 x 95 mm,
for 70 ml container, graduated, PP, translucent,
with paper label**

Page 2

Size

Sample volume	70 ml
Max. volume	70 ml
Diameter	95 mm
Diameter of opening	48,45 mm
Height of product	115 mm
Length of product	115 mm
Length excluding cap	110 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Closure material	High Density Polyethylene (HD-PE)

Purity & certification

Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	36
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	36
Piece(s) / pallet	1440
Depth of case	604 mm
Width of case	264 mm
Height of case	373 mm
Case volume	0.06 cbm
Weight of case	4.4 kg
EAN of inner box	4038917154529
EAN case	4038917036733

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT