



Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Page 1



Product description

Order number	62.547.274
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: natural, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	blue/white
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Page 2

Size

Diameter	28 mm
Length of product	114 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	300
Piece(s) / pallet	9000
Depth of case	580 mm
Width of case	330 mm
Height of case	255 mm
Case volume	0.0488 cbm
Weight of product	0.0128 kg
Weight of case	4.61 kg
EAN of inner box	4038917242561
EAN case	4038917242554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP
support and is valid without signature.



SARSTEDT