




**Product Specification**  
**Collection tip for Urine Monovette® , 100**  
**piece(s)/bag**

**Page 1**



**Product description**

Order number	10.251
Product description	Collection tip for Urine Monovette® , for Urine transfer into Monovette® Urine, 100 piece(s)/bag

**Product characteristics**

Type of collection	Aspiration
Light protection	no
Application	Urine transfer into Monovette® Urine
Suitable for	Urine Monovette®
Connection	Luer (f)

**Size**

Diameter	8 mm
Length of product	76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Collection tip for Urine Monovette® , 100**  
**piece(s)/bag**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Safety data sheet	no
Batched	yes

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	5000
Piece(s) / pallet	240000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0005 kg
Weight of case	3.31 kg
EAN of inner box	4038917080941
EAN case	4038917000116

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**