



## Product description

|                     |  |
|---------------------|--|
| Order number        | 10.252.028   |
| Product description | Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15.3 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, with plastic label, sterile, 100 piece(s)/bag |

## Product characteristics

|                       |  |
|-----------------------|--|
| Type of collection    | Aspiration                               |
| Type of preparation   | without preparation                      |
| Label/ Print          | with plastic label                       |
| Colour of print/label | transparent/yellow                       |
| Graduation            | yes                                      |
| Cap                   | Luer screw cap                           |
| Closure type          | screw cap                                |
| Light protection      | no                                       |
| Design                | Luer lock                                |
| Application           | Urine tube                               |
| Example of use        | Urine collection, transport, and storage |
| Base shape            | round base                               |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15.3 mm, 100 piece(s)/bag**

**Page 2**

|                            |           |
|----------------------------|-----------|
| Volume of work             | 10 ml     |
| Recommended centrifugation | 400 x g   |
| Centrifugation time        | 5 Min     |
| Connection                 | Luer lock |

**Size**

|                                  |         |
|----------------------------------|---------|
| Sample volume                    | 10 ml   |
| Diameter                         | 15.3 mm |
| Length including cap             | 117 mm  |
| Length excluding cap             | 102 mm  |
| Length including cap and plunger | 127 mm  |

**Material & colours**

|                         |                                   |
|-------------------------|-----------------------------------|
| Product material        | Polypropylene (PP)                |
| Colour of product       | transparent                       |
| Closure material        | High Density Polyethylene (HD-PE) |
| Colour of cap           | yellow                            |
| Tube material           | Polypropylene (PP)                |
| Piston material         | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS)                  |

**Purity & certification**

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA, ADR               |
| Product category          | no medical device   IVD |
| Purity standard           | sterile                 |
| Sterilisation             | Electron irradiation    |
| Batched                   | yes                     |

**Packaging**

|                               |            |
|-------------------------------|------------|
| Minimum order qty.            | 100        |
| Type of smallest subpackaging | bag        |
| Piece(s) / inner box          | 100        |
| Piece (s) / outer case        | 100        |
| Piece(s) / pallet             | 7000       |
| Case volume                   | 0.0229 cbm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15.3 mm, 100 piece(s)/bag**

|                   |               |
|-------------------|---------------|
| Weight of product | 0.0149 kg     |
| Weight of case    | 1.54 kg       |
| EAN of inner box  | 4038917104494 |
| EAN case          | 4038917075671 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

