



## Product Specification

**Urine Monovette® , Boric acid, 10 ml, cap green,  
(LxØ): 102 x 15 mm, 1 piece(s)/blister**

**Page 1**



### Product description

Order number	10.253.388
Product description	Urine Monovette® , Boric acid, nominal volume: 10 ml, Boric acid, Urine tube, (LxØ): 102 x 15 mm, transparent, material: PP, Luer screw cap, cap green, incl. collection tip, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix, sterile, 1 piece(s)/blister

### Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/green
Graduation	no
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Urine Monovette® , Boric acid, 10 ml, cap green,  
(LxØ): 102 x 15 mm, 1 piece(s)/blister**

**Page 2**

Volume of work	10 ml
Connection	Luer (m)

### Size

Sample volume	10 ml
Diameter	15 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	20000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500
Depth of box	351 mm
Width of box	151 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-15

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**

**Urine Monovette® , Boric acid, 10 ml, cap green,  
(LxØ): 102 x 15 mm, 1 piece(s)/blister**

Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0085 kg
Weight of case	5.89 kg
EAN of inner box	4038917950008
EAN case	4038917950015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

