




Product Specification

**V-Monovette® Urine, Boric acid, 10 ml, cap green,
(LxØ): 100 x 15 mm, 50 piece(s)/bag**

Page 1



Product description

Order number	11.2453.001
Product description	V-Monovette® Urine, Boric acid, nominal volume: 10 ml, Boric acid, vacuum system for closed urine transfer, (LxØ): 100 x 15 mm, transparent, fill mark, material: PET, membrane cap, cap green, with paper label, sterile, 50 piece(s)/bag

Product characteristics

Type of collection	Vacuum
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane cap
Closure type	push cap
Light protection	no
Application	vacuum system for closed urine transfer
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**V-Monovette® Urine, Boric acid, 10 ml, cap green,
(LxØ): 100 x 15 mm, 50 piece(s)/bag**

Page 2

Volume of work 10 ml

Size

Sample volume 10 ml
Diameter 15 mm
Length including cap 106 mm
Length excluding cap 100 mm

Material & colours

Product material Polyethylene therephthalate (PET)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green
Membrane material rubber
Tube material Polyethylene therephthalate (PET)

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet yes
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 22000
Depth of case 503 mm
Width of case 375 mm
Height of case 155 mm
Case volume 0.0292 cbm
Weight of product 0.0091 kg
Weight of case 5.02 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**V-Monovette® Urine, Boric acid, 10 ml, cap green,
(LxØ): 100 x 15 mm, 50 piece(s)/bag**

EAN of inner box 4038917263429
EAN case 4038917431903

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

