




Product Specification
Double block rack D17, Ø opening: 17 mm, 10 x 10, grey

Page 1



Product description

| | |
|---------------------|---|
| Order number | 93.852.200 |
| Product description | Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 10 x 10, grey, material: recycled PP |

Product characteristics

| | |
|------------------------|---------|
| Number of wells | 100 |
| Storage capacity | 100 |
| Format | 10 x 10 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|---------------------|--------|
| Diameter of opening | 17 mm |
| Width of product | 218 mm |
| Height of product | 45 mm |
| Length of product | 209 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Double block rack D17, Ø opening: 17 mm, 10 x 10, grey

Material & colours

| | |
|-------------------|--------------------------------------|
| Product material | recycled Polypropylene (recycled PP) |
| Colour of product | grey |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | - |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 12 |
| Piece(s) / pallet | 528 |
| Depth of case | 453 mm |
| Width of case | 433 mm |
| Height of case | 165 mm |
| Case volume | 0.0324 cbm |
| Weight of product | 0.3917 kg |
| Weight of case | 5.41 kg |
| EAN of inner box | 4038917122146 |
| EAN case | 4038917243124 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

