



Product description

Order number	93.853.134
Product description	Block Rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, blue, material: PP

Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	13 mm
Width of product	96 mm
Height of product	45 mm
Length of product	180 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	blue

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	-
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	1440
Depth of case	398 mm
Width of case	368 mm
Height of case	120 mm
Case volume	0.0176 cbm
Weight of product	0.1687 kg
Weight of case	3.16 kg
EAN of inner box	4038917122559
EAN case	4038917074056

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

