



Product Specification
**Sample tube, suitable for auto-analyzer,
transparent**

Page 1



Product description

Order number	73.641
Product description	Sample tube, suitable for auto-analyzer, material: PS, 1.5 ml, (LxØ): 24 x 14 mm, transparent, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	auto-analyzer
Base shape	conical base

Size

Max. volume	1.5 ml
Diameter	14 mm
Diameter of opening	8 mm
Height of product	24 mm
Length of product	24 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Sample tube, suitable for auto-analyzer,
transparent**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.001 kg
Weight of case	5.37 kg
EAN of inner box	4038917083072
EAN case	4038917002042

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT