



Product description

Order number	73.667
Product description	Sample tube, suitable for Abbott analyser AXSYM, material: PS, 2.5 ml, (LxØ): 53 x 12 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	Abbott analyser AXSYM
Base shape	conical base

Size

Max. volume	2.5 ml
Diameter	12 mm
Diameter of opening	10 mm
Height of product	53 mm
Length of product	53 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Sample tube, suitable for Abbott analyser AXSYM,
transparent**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2500
Piece(s) / pallet	90000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0025 kg
Weight of case	6.78 kg
EAN of inner box	4038917083232
EAN case	4038917002189

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT