



Product description

| | |
|---------------------|--|
| Order number | 73.680 |
| Product description | Mini vial, cap enclosed, material: HD-PE, 6 ml, (LxØ): 56 x 15 mm, white, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------|-----------|
| Graduation | no |
| Closure type | screw cap |
| Base shape | flat base |

Size

| | |
|---------------------|-------|
| Max. volume | 6 ml |
| Diameter | 15 mm |
| Diameter of opening | 13 mm |
| Height of product | 56 mm |
| Length of product | 56 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
| Colour of product | white |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | white |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 48000 |
| Depth of case | 570 mm |
| Width of case | 215 mm |
| Height of case | 370 mm |
| Case volume | 0.0453 cbm |
| Weight of product | 0.0033 kg |
| Weight of case | 7.4 kg |
| EAN of inner box | 4038917083348 |
| EAN case | 4038917002288 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

