



Product Specification
Multi-purpose container, 250 ml, (LxØ): 100 x 73 mm, PS, white

Page 1



Product description

Order number	75.560
Product description	Multi-purpose container, max. working volume: 250 ml, (LxØ): 100 x 73 mm, Ø opening: 73 mm, white, material: PS, 50 piece(s)/bag

Product characteristics

Graduation	no
Application	collection and storage of urine
Volume of work	250 ml

Size

Sample volume	250 ml
Diameter	73 mm
Diameter of opening	73 mm
Height of product	100 mm
Length of product	100 mm
Length excluding cap	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Multi-purpose container, 250 ml, (LxØ): 100 x 73 mm, PS, white

Material & colours

Product material	Polystyrene (PS)
Colour of product	white

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	20000
Depth of case	598 mm
Width of case	398 mm
Height of case	165 mm
Case volume	0.0393 cbm
Weight of product	0.0076 kg
Weight of case	4.3 kg
EAN of inner box	4038917339902
EAN case	4038917339896

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

