



**Product Specification**  
**Medicine cup, 50 ml, (LxØ): 54 x 55 mm,**  
**graduated, PP, transparent**

**Page 1**



**Product description**

Order number	75.9924.522
Product description	Medicine cup, max. working volume: 50 ml, (LxØ): 54 x 55 mm, Ø opening: 55.4 mm, transparent, graduated, material: PP

**Product characteristics**

Graduation	yes
Application	Administration of medicines
Example of use	Oral administration of medications in liquid or tablet form
Volume of work	50 ml
Autoclavable (max. °C)	121 °C

**Size**

Sample volume	50 ml
Max. volume	50 ml
Diameter	55 mm
Diameter of opening	55,4 mm
Height of product	54 mm
Length of product	54 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Medicine cup, 50 ml, (LxØ): 54 x 55 mm,**  
**graduated, PP, transparent**

**Page 2**

Length excluding cap 55.4 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category Medical device  
CE certified CE  
Latex-free yes  
Pyrogen-free no  
Safety data sheet no

**Packaging**

Minimum order qty. 500  
Piece(s) / inner box 500  
Piece (s) / outer case 500  
Piece(s) / pallet 80000  
Depth of case 271 mm  
Width of case 224 mm  
Height of case 158 mm  
Case volume 0.0099 cbm  
Weight of case 3 kg  
EAN of inner box 4038917352949  
EAN case 4038917352932

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**