



Product Specification
Tube segment opener, mounted on tube
75 x 12 mm, HD-PE, white, tube: PS

Page 1



Product description

Order number	95.1000.476
Product description	Tube segment opener, mounted on tube 75 x 12 mm, material: HD-PE, white, tube: PS, 50 piece(s)/StackPack

Product characteristics

Design	mounted on tube 75 x 12 mm
Application	Blood bank
Base shape	round base

Size

Diameter	12 mm
Length of product	75 mm

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	white
Tube material	Polystyrene (PS)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification
Tube segment opener, mounted on tube
75 x 12 mm, HD-PE, white, tube: PS

Page 2

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	400
Type of smallest subpackaging	StackPack
Piece(s) / inner box	50
Piece (s) / outer case	400
Piece(s) / pallet	57600
Depth of box	146 mm
Width of box	96 mm
Height of box	64 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0045 kg
Weight of case	2.22 kg
EAN of inner box	4038917124614
EAN case	4038917047975

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT