



Product Specification

Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,
conical false bottom, flat tube bottom, PP, cap
assembled, 100 piece(s)/bag

Page 1



Product description

Order number	62.611
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, transparent, material: PP, with print, label/print: white, with graduations, cap assembled, 100 piece(s)/bag, 1,000 piece(s)/carton

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Base shape	conical false bottom, flat tube bottom
Volume of work	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,
conical false bottom, flat tube bottom, PP, cap
assembled, 100 piece(s)/bag**

Page 2

Size

Diameter	15.3 mm
Length of product	92 mm
Length excluding cap	92 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.005 kg
Weight of case	5.69 kg
EAN of inner box	4038917082877
EAN case	4038917001892

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT