



### **Product description**

Order number	65.1379.006
Product description	Cap, blue, suitable for tubes Ø 12-17 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard cap
Suitable for	tubes Ø 12-17 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	blue

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cap, blue, suitable for tubes Ø 12-17 mm**

**Page 2**

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	240000
Depth of case	413 mm
Width of case	257 mm
Height of case	464 mm
Case volume	0.0492 cbm
Weight of product	0.0004 kg
Weight of case	5.36 kg
EAN of inner box	4038917128018
EAN case	4038917021548

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**