



Product description

Order number	62.610.018
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 92 x 15.3 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, with print, label/print: green, "CSF/Liquor", with graduations, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with print
Colour of print/label	green
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	10 ml
Centrifugation max (RCF)	3000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 10 ml, (LxØ): 92 x 15.3 mm, PP,
with print**

Page 2

Size

Diameter	15.3 mm
Length of product	92 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500
Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0052 kg
Weight of case	4.23 kg
EAN of inner box	4038917233132
EAN case	4038917233101

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

