



Product description

Order number 73.834.001
Product description Steel ball for Amelung KC coagulometer (macro), grey

Product characteristics

Suitable for Amelung KC coagulometer (macro)

Material & colours

Colour of product grey

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Piece(s) / inner box 1000
Piece (s) / outer case 50000
Piece(s) / pallet 5400000
Depth of case 308 mm
Width of case 223 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Steel ball for Amelung KC coagulometer (macro)

Page 2

Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0001 kg
Weight of case	7.15 kg
EAN of inner box	4038917121682
EAN case	4038917017930

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT