



### **Product description**

Order number	86.1996
Product description	ESR pipette, with fitted piston insert

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of box	233 mm
Width of box	164 mm
Height of box	76 mm
Depth of case	388 mm
Width of case	338 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**ESR pipette**

**Page 2**

Height of case	250 mm
Case volume	0.0328 cbm
Weight of product	0.0076 kg
Weight of case	8.61 kg
EAN of inner box	4038917088138
EAN case	4038917006231

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**