



Product description

Order number	90.1060.062
Product description	ESR rack, suitable for S-Monovette® , without back panel

Product characteristics

Design	without back panel
Suitable for	S-Monovette®

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	1560
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
ESR rack, suitable for S-Monovette®

Page 2

Case volume	0.018 cbm
Weight of product	1.167 kg
Weight of case	23.8 kg
EAN of inner box	4038917124935
EAN case	4038917019545

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT