




Product Specification

ESR rack, suitable for Microvette[®] , 10 measuring stations

Page 1



Product description

Order number	90.1091
Product description	ESR rack, suitable for Microvette [®] , 10 measuring stations

Product characteristics

Suitable for	Microvette [®]
Number of positions	10

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	8
Piece(s) / pallet	160
Depth of case	263 mm
Width of case	98 mm
Height of case	335 mm
Case volume	0.0086 cbm
Weight of product	0.2505 kg
Weight of case	5.1 kg
EAN of inner box	4038917085311
EAN case	4038917004015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

