



Product Specification
**Salivette® , with cotton swab, cap: white, with
paper label**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 51.1534 |
| Product description | Salivette® , with cotton swab, cap: white, with paper label, label/print: white/black, for Saliva examination, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|--------------------|
| Label/ Print | with paper label |
| Colour of print/label | white/black |
| Closure type | push cap |
| Design | with cotton swab |
| Application | Saliva examination |
| Recommended centrifugation | 1000 x g |
| Centrifugation time | 2 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|----------------------|---------|
| Diameter | 16.8 mm |
| Length including cap | 97 mm |
| Length excluding cap | 84 mm |

Material & colours

| | |
|------------------|----------------------------------|
| Product material | Polypropylene (PP) |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | white |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 28000 |
| Depth of case | 598 mm |
| Width of case | 398 mm |
| Height of case | 125 mm |
| Case volume | 0.0298 cbm |
| Weight of product | 0.007 kg |
| Weight of case | 4.09 kg |
| EAN of inner box | 4038917086905 |
| EAN case | 4038917005432 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

