



**Product Specification**  
**Polystyrene container, suitable as outer packaging for cool transport**

**Page 1**



**Product description**

Order number	95.1011
Product description	Polystyrene container, suitable as outer packaging for cool transport

**Product characteristics**

Design	outer packaging for cool transport
Absorbent liner	no

**Size**

Width of product	71 mm
Height of product	168 mm
Length of product	140 mm

**Material & colours**

Product material	styrofoam
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Polystyrene container, suitable as outer  
packaging for cool transport**

**Page 2**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      1  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      1  
Piece (s) / outer case                      36  
Piece(s) / pallet                              720  
Depth of case                              584 mm  
Width of case                                387 mm  
Height of case                               360 mm  
Case volume                                0.0814 cbm  
Weight of product                        0.0371 kg  
Weight of case                              2.25 kg  
EAN of inner box                        4038917084987  
EAN case                                    4038917362399

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**