



Product description

Order number	95.1710
Product description	Foam film bag, (LxW): 185 x 230 mm, material: PE

Size

Width of product	230 mm
Length of product	185 mm

Material & colours

Product material	Polyethylene (PE)
------------------	-------------------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	200
Piece(s) / inner box	200
Piece (s) / outer case	200
Piece(s) / pallet	6400
Depth of case	458 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Foam film bag, (LxW): 185 x 230 mm, PE

Page 2

Width of case	408 mm
Height of case	205 mm
Case volume	0.0383 cbm
Weight of product	0.0018 kg
Weight of case	0.91 kg
EAN of inner box	4038917087490
EAN case	4038917079648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT