



**Product description**

Order number	95.1717.011
Product description	Transport case B 13, flat design for block rack, (LxWxH): 395 x 395 x 300 mm

**Product characteristics**

Design	flat design for block rack
--------	----------------------------

**Size**

Width of product	395 mm
Height of product	300 mm
Length of product	395 mm

**Material & colours**

Product material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	30
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	2.763 kg
Weight of case	3.3 kg
EAN of inner box	4038917133197
EAN case	4038917069250

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

