



Product Specification

Transport case I 100, internal case with attachment

Page 1



Product description

Order number	95.1720
Product description	Transport case I 100, internal case with attachment, (LxWxH): 360 x 310 x 275 mm

Product characteristics

Design	internal case with attachment
Shape of inner tube	square

Size

Width of product	310 mm
Height of product	275 mm
Length of product	360 mm

Material & colours

Product material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Satisfies the requirement	TRBA 100
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	30
Depth of case	398 mm
Width of case	358 mm
Height of case	325 mm
Case volume	0.0463 cbm
Weight of product	1.613 kg
Weight of case	2.05 kg
EAN of inner box	4038917087551
EAN case	4038917076210

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

