



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 95.1725   |
| Product description | Battery bag, without batteries, (LxW): 205 x 115 mm, material: PL, for standard cooling batteries, 12 piece(s)/case |

### **Product characteristics**

|              |                            |
|--------------|----------------------------|
| Design       | without batteries          |
| Suitable for | standard cooling batteries |

### **Size**

|                   |        |
|-------------------|--------|
| Width of product  | 115 mm |
| Length of product | 205 mm |

### **Material & colours**

|                  |                |
|------------------|----------------|
| Product material | polyester (PL) |
|------------------|----------------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 12  
Type of smallest case  
subpackaging  
Piece(s) / inner box 12  
Piece (s) / outer case 12  
Piece(s) / pallet 1980  
Depth of case 338 mm  
Width of case 228 mm  
Height of case 120 mm  
Case volume 0.0092 cbm  
Weight of product 0.4419 kg  
Weight of case 5.5 kg  
EAN of inner box 4038917087568  
EAN case 4038917068376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

