



Product Specification
Screw cap, natural, suitable for Mailing container
126 x 30 mm

Page 1



Product description

| | |
|---------------------|---|
| Order number | 65.679 |
| Product description | Screw cap, natural, suitable for Mailing container 126 x 30 mm, 50 piece(s)/bag |

Product characteristics

| | |
|--------------|-------------------------------|
| Application | Mailing of samples |
| Suitable for | Mailing container 126 x 30 mm |

Size

| | |
|-------------------|-------|
| Length of product | 20 mm |
|-------------------|-------|

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | natural |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap, natural, suitable for Mailing container
126 x 30 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 250
Type of smallest subpackaging bag
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 41250
Depth of case 338 mm
Width of case 228 mm
Height of case 120 mm
Case volume 0.0092 cbm
Weight of product 0.0035 kg
Weight of case 1.1 kg
EAN of inner box 4038917083331
EAN case 4038917002271

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

