



Product Specification

Mailing container, length: 179 mm, Ø opening: 28 mm, cap enclosed

Page 1



Product description

| | |
|---------------------|--|
| Order number | 78.574 |
| Product description | Mailing container, (LxW): 179 x 35 mm, Ø opening: 28 mm, cap enclosed, cap: natural, 50 piece(s)/bag |

Product characteristics

| | |
|------------------|--|
| Cap | round cap |
| Closure type | screw cap |
| Application | Mailing of samples |
| Example of use | Secondary packaging according to P 650 |
| Absorbent liner | no |
| Shape of opening | round |
| Liquid-tight | yes |

Size

| | |
|---------------------|--------|
| Diameter of opening | 28 mm |
| Width of product | 35 mm |
| Length of product | 179 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
| Colour of cap | natural |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650, ADR |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 250 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 250 |
| Piece(s) / pallet | 3000 |
| Depth of case | 563 mm |
| Width of case | 343 mm |
| Height of case | 500 mm |
| Case volume | 0.0966 cbm |
| Weight of product | 0.0201 kg |
| Weight of case | 5.94 kg |
| EAN of inner box | 4038917159449 |
| EAN case | 4038917159432 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

