



Product Specification
**Absorbent liner, suitable for Mailing container 114
x 44 mm, (LxW): 135 x 95 mm**

Page 1



Product description

Order number 78.669
Product description Absorbent liner, suitable for Mailing container 114 x 44 mm,
(LxW): 135 x 95 mm, 100 piece(s)/case

Product characteristics

Application Mailing of samples, Multi-Safe disposel boxes
Suitable for Mailing container 114 x 44 mm

Size

Width of product 95 mm
Length of product 135 mm

Material & colours

Product material cellulose

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Absorbent liner, suitable for Mailing container 114
x 44 mm, (LxW): 135 x 95 mm

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	385 mm
Width of case	385 mm
Height of case	160 mm
Case volume	0.0237 cbm
Weight of product	0.0026 kg
Weight of case	3.04 kg
EAN of inner box	4038917083249
EAN case	4038917002196

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

